






# Corning® HYPERFlask® Cell Culture Vessel jetPEI™ Fast Transfection Protocol

## Protocol

This is an optimized rapid version of the standard jetPEI High Yielding Performance *Flask* (HYPER*Flask*) cell culture vessel transfection protocol (CLS-AN-110). It was developed to save time and cost by transfecting and seeding the HYPER*Flask* vessel in one step.

### Introduction

One of the most useful tools in cell biology research is transfection, the introduction of foreign DNA into eukaryotic cells. In much of today's research, there is a growing need for the effective transfection of large quantities of cells. Polyplus-transfection's jetPEI transfection reagent is a highly efficient, low toxicity, water-soluble polymer that can be used in the presence of serum in culture media. Therefore, there is no need to change the culture medium before or after transfecting cells, making this method ideally suited for use with the Corning HYPER*Flask* cell culture vessel. This protocol was optimized using HeLa cells but has been successfully applied to a variety of cell types including Chinese hamster ovary (CHO) cells. This protocol is intended as a starting point that can be optimized by the end user for their cell lines.



## Cell Culture Suspension Set Up

Steps were modified from Corning HYPERFlask® Cell Culture Vessel, jetPEI™ Transfection Protocol and Polyplus-transfection's jetPEI™ transfection protocol.

The procedure described below is for plating cells into a HYPERFlask cell culture vessel (Corning Cat. No. 10010 or 10024) and multiple wells of a 24 well plate (Corning Cat. No. 3524). The 24 well plate will serve as a control for overall transfection efficiency as well as the transfection efficiency of the large scale precipitate made for the HYPERFlask vessel. Should you choose to use a different size control well, scale your changes in reagents based on an equivalent mL/cm<sup>2</sup> (Table 1).

*Helpful Hint:* Use early passage cultures (5 to 20 passages) at 80 to 90% confluence.

1. Harvest flask using standard harvest techniques, centrifuge for 10 min at 1,000 rpm and re-suspend cells in 50 mL of fresh growth medium.

*Helpful Hint:* If necessary, pass cell suspension through cell strainer to achieve a single cell suspension.

2. Prepare a cell suspension at  $1.82 \times 10^5$  cells/mL in a final volume of 424 mL for each HYPERFlask vessel to be tested using fresh growth medium.
3. Prepare a cell suspension at  $1.82 \times 10^5$  cells/mL in a final volume of 550  $\mu$ L for each control or mock well to be tested.

*Note:* All mock or control replicates can be made as one cocktail and split over each well.

## JetPEI/DNA Complex

Steps have been modified from Polyplus' jetPEI transfection protocol using jetPEI transfection kit (Part No. 101-40N) optimized for 1  $\mu$ g DNA and a jetPEI N:P ratio of 5.

1. DNA solution A. Prepare in a container or tube that can hold 2X the final volume.

Solution A*	For One 24 Well Mock (0.650 mL/Well)	For One 24 Well Control (0.650 mL/Well)	For One HYPERFlask Vessel (560 mL/Flask)
DNA	–	0.5 $\mu$ g/cm <sup>2</sup> (1 $\mu$ g)	0.5 $\mu$ g/cm <sup>2</sup> (860 $\mu$ g)
150 mM NaCl	50 $\mu$ L	To 50 $\mu$ L	To 43.12 mL
<b>Final Volume</b>	<b>50 <math>\mu</math>L</b>	<b>50 <math>\mu</math>L</b>	<b>43.1 mL</b>

\*Prepare in a container/tube that can hold 2X the final volume.

2. JetPEI™ Solution B

Solution B**	For One 24 Well Mock (0.650 mL/Well)	For One 24 Well Control (0.650 mL/Well)	For One HYPERFlask Vessel (560 mL/Flask)
JetPEI Reagent	2 $\mu$ L	2 $\mu$ L	1.72 mL
150 mM NaCl	48 $\mu$ L	48 $\mu$ L	41.4 mL
<b>Final Volume</b>	<b>50 <math>\mu</math>L</b>	<b>50 <math>\mu</math>L</b>	<b>43.1 mL</b>

\*\*Optimized for an N:P ratio of 5.

3. Rapidly, while gently mixing add jetPEI solution B into DNA solution A; mix well by vortexing.

**Important Note: Do not add in reverse order.**

*Note:* All mock or control replicates can be made as one cocktail and split over each well.

Final Volume	For One 24 Well	For One HYPERFlask Vessel
JetPEI/DNA Complex	100 $\mu$ L	86.2 mL

4. Incubate at room temperature for 15 min. Solution may appear cloudy.

## Transfection

For handling of the HYPERFlask® cell culture vessel refer to the HYPERFlask Cell Culture Vessel Instructions for Use.

### HYPERFlask Cell Culture Vessel

1. Slowly, while mixing, add 86.2 mL jetPEI™/DNA complex to the cell suspension.

*Note:* Minimize foaming of medium when mixing.

2. Seed 650 µL/well of cell/DNA complex solution to 24 well control plate to serve as the large scale precipitate controls.

*Note:* Up to 3 wells can be tested for performance of the large scale complex in 24 well plate without interfering with the efficiency of transfection of the HYPERFlask vessel.

3. Gently pour remaining solution into HYPERFlask vessel, remove all trapped air and recap securely

*Helpful Hint:* If necessary, use fresh growth medium to bring liquid volume in the HYPERFlask vessel to the first thread.

### Control Wells

4. Slowly, while gently mixing add the jetPEI/DNA complex to cell suspension.
5. Seed 650 µL/well of the correct solution into corresponding control wells.
6. Incubate all cultures at 5% CO<sub>2</sub>, 37°C for 48 hours.
7. Process transfected cells as necessary.

Please visit the Corning Life Sciences web site to view a video presentation that describes the proper handling of the HYPERFlask cell culture vessel.

For additional product or technical information, please e-mail us at [CLStechserv@corning.com](mailto:CLStechserv@corning.com), visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences), or call 1.800.492.1110. Outside the United States, please call 1.978.442.2200.

# CORNING

## **Corning Incorporated** *Life Sciences*

Tower 2, 4th Floor  
900 Chelmsford St.  
Lowell, MA 01851  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

## **Worldwide Support Offices**

### **ASIA / PACIFIC**

**Australia/New Zealand**  
t 65 6733-6511  
f 65 6735-2913

**China**  
t 86-21-5467-4666  
f 86-21-5407-5899

**India**  
t 91-124-235 7850  
f 91-124-401 0207

**Japan**  
t 81 3-3586 1996  
f 81 3-3586 1291

**Korea**  
t 82 2-796-9500  
f 82 2-796-9300

**Singapore**  
t 65 6733-6511  
f 65 6861-2913

**Taiwan**  
t 886 2-2716-0338  
f 886 2-2716-0339

## **EUROPE**

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 31 20 655 79 28  
f 31 20 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

## **All Other European Countries**

t 31 (0) 20 659 60 51  
f 31 (0) 20 659 76 73

## **LATIN AMERICA**

**Brasil**  
t (55-11) 3089-7419  
f (55-11) 3167-0700

**Mexico**  
t (52-81) 8158-8400  
f (52-81) 8313-8589