

Robotic Arm Stackers and Plate Movement - Compatible Corning Products		96 Well Plates								384 Well Plates					1536 Well Plates				Polypropylene Plates									
Company Name	Product Name	96 well flat bottom	96 well round bottom	96 well V bottom	96 well Easy-Wash bottom	96 well half area flat bottom	96 well clear flat bottom	96 well half area clear flat bottom	96 well special optics	384 well flat bottom	384 well low volume round bottom	384 well w/clear bottom	384 well optical imaging	384 well low volume clear flat bottom	1536 10uL	1536 well black clear bottom	1536 well flat bottom	1536 2uL "wafer"	96 well flat bottom PP	96 well round bottom PP	96 well V bottom PP	Clear Pro Round PP	96 .5 mL V Block PP	96 well 1 mL round bottom block PP	96 well 2ml Block PP	384 well Plate PP	384 well low volume round bottom PP	384 well round block PP
Beckman Coulter	ORCA Robot®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓				
Bio-Tek® Instruments, Inc.	BioStack™	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓						✓	✓	✓	✓	✓			✓	✓	
Caliper Life Sciences, Inc.	Allegro™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓	✓	✓	✓				
Caliper Life Sciences, Inc.	Presto® AutoStack	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Caliper Life Sciences, Inc.	Twister® II	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
CyBio	CyBio Lifter	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓	✓	
CyBio	CyBio Rotating arm	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓	✓	
CyBio	CyBio Stackers	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓	✓
Hudson Control Group	PlateCrane™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Hudson Control Group	PlateCrane™ XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Hudson Control Group	PlateCrane™ XT	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Hudson Control Group	PlateSilo™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Kendro	Teleshake	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓	✓	
LabCyte	MiniStack	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓			✓	✓	
PerkinElmer Life & Analytical Sciences	MiniTrak™	✓	✓	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					✓
PerkinElmer Life & Analytical Sciences	PlateStak™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
PerkinElmer Life & Analytical Sciences	PlateTrak™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓			
Tecan	Mol Bank	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Tecan	Tecan Carousel	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Tecan	TRAC™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
TekCel	PlateServer Ultra High Capacity Storage TekBench	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓	✓	

