



AUTOMATION - Plate Sealers

Company Name	Product Name	Image	Plate Height			Bottom Side	Top Side	Heat Seal	Pressure Sensitive	Integratable
			Small	Medium	Deep					
AB Gene®	ALPS Plate Sealer		Yes	Yes	Yes	No	Yes	Yes	Yes	Yes
AB Gene®	Thermosealer		Yes	Yes	Yes	No	Yes	Yes	No	No
Beckman Coulter	SAGIAN™ Sealer		Yes	Yes	Yes	No	Yes	No	Yes	Yes
Tomtec	Quadra® Seal		Yes	Yes	Yes	Yes	Yes	No	Yes	Yes
Velocity 11	PlateLoc® Thermal Plate Sealer		Yes	Yes	Yes	No	Yes	No	Yes	Yes