

Microplate Dispensers - Compatible Corning Products		96 Well Plates								384 Well Plates					1536 Well Plates				Polypropylene Plates									
Company Name	Product Name	96 well flat bottom	96 well round bottom	96 well V bottom	96 well Easy-Wash bottom	96 well half area flat bottom	96 well clear flat bottom	96 well half area clear flat bottom	96 well Special Optics	384 well flat bottom	384 well low volume round bottom	384 well w/clear bottom	384 well Optical Imaging	384 well low volume clear flat bottom	1536 10uL	1536 well black clear bottom	1536 well flat bottom	1536 2uL "wafer"	96 well flat bottom PP	96 well round bottom PP	96 well V bottom PP	Clear Pro Round PP	96 .5 mL V Block	96 well 1 mL round bottom block	96 well 2mL Block	384 well plate	384 well black low volume round bottom PP	384 well round block
Beckman Coulter, Inc.	BioRAPTR FDR™	X	X	X	X	X	X	X	X	X		X	X	X		X	X											
Artel	Artel MVS (Multichannel Verification System)	X	X	X	X	X	X	X	X	X	X	X	X	X					X	X	X	X	X	X	X	X	X	X
BioMicroLab, Inc	VolumeCheck	X	X	X	X	X	X	X	X										X	X	X	X	X	X	X			
Bio-Tek® Instruments, Inc	µFill™	X	X	X	X	X	X	X	X	X		X	X	X														
Bio-Tek® Instruments, Inc	Precision™	X	X	X	X	X	X	X	X	X	X	X	X	X					X	X	X	X	X			X	X	X
Bio-Tek® Instruments, Inc	Precision™ XS	X	X	X	X	X	X	X	X	X	X	X	X	X					X	X	X	X	X			X	X	X
Brandel	PXR-96-MS	X	X	X	X	X	X	X	X										X	X	X	X	X	X	X			
Deerac Fluidics™	Equator™ HTS	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	X
Genetix	Qfill 2	X	X	X	X	X	X	X	X	X	X	X	X	X														
Innovadyne Technologies, Inc.	Nanodrop™ I	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	X
Innovadyne Technologies, Inc.	Nanodrop™ II	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	X
Innovadyne Technologies, Inc.	Screenmaker96+8™	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	X
Labcyte	MiniGene Line Dispenser	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X											
Perkin Elmer Life & Analytical Science	FlexDrop™ III	X	X	X	X	X	X	X	X	X		X	X	X		X	X											
Thermo Electron Corporation	Multidrop™ 384	X	X	X	X	X	X	X	X	X	X	X	X	X					X	X	X	X						
Thermo Electron Corporation	Multidrop™ 96	X	X	X	X	X	X	X	X										X	X	X	X						
Thermo Electron Corporation	Multidrop™ Deep Well	X	X	X	X	X	X	X	X										X	X	X	X	X	X	X			