



LIQUID HANDLING - Acoustic Dispensers										
Company Name	Product Name	Image	Microplate Formats			Dynamic Range	Dispensing Precision	Accuracy	Resolution	Throughput
			96	384	1536					
Labcyte	Echo®520 Echo®550 Echo®555		Yes	Yes	Yes	2.5nL - >5mL	<8%CV	<10% Deviation from target	2.5nL	<1.5Min./384 single drop transfer
EDC Biosystems	ATS-100		Yes	Yes	Yes	0.5 - 10µL+	<10% CV	<10% Deviation from target	1-10 nL	<6 Min./1536 single drop transfer