



## PLATE READERS - Luminescence

Company Name	Product Name	Image	Microplate Formats			Luminescence Mode		Top Read	Bottom Read	Wavelength Range	Injectors	Incubation	Shaking	Integratable
			96	384	1536	Flash	Glow							
Applied Biosystems	TR717™		Yes	Yes	No	Yes	Yes	Yes	No	380-630 nm	0-2	Yes	No	Yes
Berthold Technologies GmbH & Co. KG	Centro LB 960		Yes	Yes	No	Yes	Yes	Yes	No	380-630 nm	0-3	Yes	Yes	Yes
Berthold Technologies GmbH & Co. KG	MPL1		Yes	No	No	No	Yes	Yes	No	300-600 nm	No	No	No	No
Berthold Technologies GmbH & Co. KG	MPL2		Yes	No	No	No	Yes	Yes	No	300-600 nm	No	No	No	No
Berthold Technologies GmbH & Co. KG	ORION L		Yes	Yes	No	Yes	Yes	Yes	No	300-600 nm	0-2	No	Yes	No
Berthold Technologies GmbH & Co. KG	ORION II		Yes	Yes	No	Yes	Yes	Yes	No	300-600 nm	0-4	Yes	Yes	Yes

## PLATE READERS - Luminescence

Company Name	Product Name	Image	Microplate Formats			Luminescence Mode		Top Read	Bottom Read	Wavelength Range	Injectors	Incubation	Shaking	Integratable
			96	384	1536	Flash	Glow							
Bio-Tek® Instruments	Clarity™		Yes	Yes	No	Yes	Yes	Yes	No	300-600nm	Yes/ Dual	Yes	Yes	Yes
BMG LABTECH	LUMIstar OMEGA		Yes	Yes	Yes	Yes	Yes	Yes	Yes	200-900 nm	Yes	Yes	Yes	Yes
Harta Instruments Inc.	MicroLumi XS		Yes	No	No	Yes	Yes	Yes	No	280-600 nm	up to 2	No	No	Yes
Thermo Electron Corporation	Luminoskan Ascent		Yes	Yes	No	Yes	Yes	Yes	Yes	185-680 nm	0-3	Yes	Yes	Yes
Turner BioSystems, Inc.	Veritas™		Yes	No	No	Yes	Yes	Yes	No	350-650 nm	0-2	No	No	No

**PLATE READERS - Luminescence**

Company Name	Product Name	Image	Microplate Formats			Luminescence Mode		Top Read	Bottom Read	Wavelength Range	Injectors	Incubation	Shaking	Integratable
			96	384	1536	Flash	Glow							
Turner BioSystems, Inc.	GloRunner™		Yes	No	No	No	Yes	Yes	No	300-610 nm	No	No	No	No