

### PLATE READERS - Scintillation Counters

Company Name	Product Name	Image	Microplate Formats			Scintillation Counting	Luminescence Mode		Injectors	# of Detectors	Top and/or Bottom Read	Integratable
			96	384	1536		Flash	Flow				
PerkinElmer Life & Analytical Sciences	MicroBeta® JET		Yes	Yes	No	Yes	Yes	Yes	Yes	Up to 6	Top or bottom read for Lumi. Coincidence counting for SPA & LSC.	Yes
PerkinElmer Life & Analytical Sciences	MicroBeta® TriLux		Yes	Yes	No	Yes	No	Yes	No	Up to 12	Top or bottom read for Lumi. Coincidence counting for SPA & LSC.	Yes
PerkinElmer Life & Analytical Sciences	TopCount® NXT		Yes	No	No	Yes	No	Yes	No	Up to 12	Top read	Yes
PerkinElmer Life & Analytical Sciences	TopCount® NXT-HTS		Yes	Yes	No	Yes	No	Yes	No	Up to 12	Top read	Yes