

Absorbance - Compatible Corning Products		96 Well Plates								384 Well Plates					1536 Well Plates				Polypropylene Plates									
Company Name	Product Name	96 well flat bottom	96 well round bottom	96 well V bottom	96 well Easy-Wash bottom	96 well half area flat bottom	96 well clear flat bottom	96 well half area clear flat bottom	96 well special optics	384 well flat bottom	384 well low volume round bottom	384 well w/clear bottom	384 well optical imaging	384 well low volume clear flat bottom	1536 10uL	1536 well black clear bottom	1536 well flat bottom	1536 2uL "wafer"	96 well flat bottom PP	96 well round bottom PP	96 well V bottom PP	Clear Pro Round PP	96 .5 mL V Block PP	96 well 1 mL round bottom block PP	96 well 2ml Block PP	384 well Plate PP	384 well low volume round bottom PP	384 well round block
Anthos Labtec	2010	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Anthos Labtec	2020	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Anthos Labtec	Zenyth 340	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
Anthos Labtec	Zenyth 200	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓						
Bio-Rad	Model 550	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Bio-Rad	Ultramar™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						
Bio-Rad®	Benchmark™	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Bio-Tek® Instruments	PowerWave™	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Bio-Tek® Instruments	uQuant™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓						
Bio-Tek® Instruments	ELx800™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓					✓	✓	✓	✓						
Bio-Tek® Instruments	ELx808™	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
BMG LABTECH	SPECTROstar	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓				✓	✓	
Molecular Devices Corporation	SpectraMax® Plus 384	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓														
Molecular Devices Corporation	SpectraMax® 340 PC 384	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓														
Molecular Devices Corporation	SpectraMax® 190	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Molecular Devices Corporation	SpectraMax® 340PC	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Molecular Devices Corporation	SpectraMax® Plus	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Molecular Devices Corporation	EMax®	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						
Molecular Devices Corporation	VMax®	✓	✓	✓	✓	✓	✓	✓	✓										✓	✓	✓	✓						

